

# EquaGel Adjustable Protector

---

<https://vital.net.au/products/pressure-care-skin-protection/sitting-pressure-care/equagel-adjustable-protector>



The EquaGel Adjustable Protector Cushion is an alternative cushion to any other on the market. Its unique gel design offers maximum comfort by distributing weight evenly when seated. It has been designed to effectively lower and minimize the pressure on the buttocks reducing the chances of pressure sores. It is a great product for people who have very limited movement while seated and are prone to pressure ulcers.

This pressure cushion is made by two layers of dry-polymer gel; the top layer comprises a thin wall of gel, while the lower is made from a thicker gel layer, designed to maintain posture support. The Adjustable Protector cushion is 2.5" in height and is contoured in the pelvic area to more evenly distribute weight. The Adjustable cushion has two gel inserts to customise the cushion to the user. The inserts are 1/4" and 1/2" and can be used separately or together to adjust the height under the pressure areas from 1 3/4" up to 2", 2 1/4" or 2 1/2" deep.

The EquaGel Adjustable Protector Cushion has a waterproof cover and velcro strip on the bottom for wheelchair users.

**Product Code:** 8-400027-1

**Price:** \$679.40

## Features:

The EquaGel Adjustable Protector Cushion is an alternative cushion to any other on the market. Its unique gel design offers maximum comfort by distributing weight evenly when seated. It has been designed to effectively lower and minimize the pressure on the buttocks reducing the chances of pressure sores. It is a great product for people who have very limited movement while seated and are prone to pressure ulcers.

This pressure cushion is made by two layers of dry-polymer gel; the top layer comprises a thin wall of gel, while the lower is made from a thicker gel layer, designed to maintain posture support. The Adjustable Protector cushion is 2.5" in height and is contoured in the pelvic area to more evenly distribute weight. The Adjustable cushion has two gel inserts to customise the cushion to the user. The inserts are 1/4" and 1/2" and can be used separately or together to adjust the height under the pressure areas from 1 3/4" up to 2", 2 1/4"

or 2 1/2? deep.

The EquaGel Adjustable Protector Cushion has a waterproof cover and velcro strip on the bottom for wheelchair users.

### **Specifications:**

- The EquaGel Adjustable Protector Cushion has a waterproof cover and velcro strip on the bottom for wheelchair users.
- Multiple sizes

### **Shipping Information:**

**Delivery is available throughout Australia** for all online purchases. **Free In-Store Pick Up** is available from our Forster, Taree and Port Macquarie showrooms.

Delivery fees for small items start at \$12 and vary depending on the product size and delivery destination. Fees are calculated in the Review & Payment section prior to checkout.

Delivery timeframe is subject to product availability.

**If you require your item(s) urgently** please contact Vital Living on 02 6581 2400 or via email [customerservice@vital.net.au](mailto:customerservice@vital.net.au) prior to completing your online purchase to **check stock availability** and estimated delivery timeframe.

Free delivery is available for certain purchases over \$250 for online orders only. **Conditions apply\*** depending on delivery destination and you will be notified of any delivery fees in the Review & Payment section prior to checkout.

[Click here for more detailed shipping information](#)